
Elis 1C

Model: R-A033

RENPHO is committed to our environment. To reduce paper usage, please scan the QR Code or enter the following URL to download the full version of the electronic User Manual:

<https://renpho.com/pages/product-manual-body-fat-scale>

Contents

English..... EN-01~EN-10

Deutsch..... DE-11~DE-20

Français..... FR-21~FR-29

Italiano..... IT-30~IT-38

Español..... ES-39~ES-47

Thank you for choosing RENPHO Elis 1C. This scale offers you detailed insights into your 13 health metrics to help achieve a healthier lifestyle. We wish you a pleasant experience and success on your health and fitness journey!

IMPORTANT SAFETY INSTRUCTIONS

READ ALL INSTRUCTIONS CAREFULLY BEFORE USE.

1. The scale is for household use only and not for medical or commercial purposes.
2. Children under 8 or individuals who are pregnant, have osteoporosis, a pacemaker, artificial joints, portable ECG or other metal implants should use the scale for weighing purposes only and make sure to wear shoes to prevent any electric current when stepping onto the scale.
3. All body composition data provided is for reference and should not be used as a substitute for advice from a licensed healthcare professional.
4. Do not overload the scale. (Maximum Capacity: 440lb / 200kg)
5. Inspect the scale for any damage before each use. Do not use a damaged scale to avoid inaccuracy or injury.
6. In case of any malfunction, check the battery power first. If necessary, replace the batteries.
7. Keep the batteries out of reach of children. In case of ingestion of batteries, seek medical assistance promptly.
8. To ensure accuracy, use the scale on a hard, flat and dry surface.
9. To avoid interference or danger, do not use the scale in an environment near equipment with strong electromagnetic fields.
10. To prevent injury, do not stand on the edges of the scale, jump on or stomp it.
11. The scale is a delicate and precise measuring device. Handle with care to avoid breakage.
12. The scale is not waterproof. Keep it dry and never submerge it in water. To clean the surface, use a damp cloth or glass cleaner. Do not use soap or chemical cleaning agents.
13. Do not expose the scale to direct sunlight or use it in high temperatures or extremely humid areas. Store the scale in a cool, dry place.

14. Do not attempt to disassemble the scale for any reason, as this may cause damage or void the warranty.
15. Exhausted batteries should be immediately removed from the scale and properly disposed of. Do not dispose of batteries in fire. When discharged batteries are kept in the scale for a long time, electrolyte leakage can occur causing damage to the scale and/or personal injury.
16. Do not heat, short-circuit, force discharge, encapsulate and/or modify the batteries. Do not deform or dismantle the batteries.
17. Do not weld or solder directly to batteries.
18. Store unused batteries in their original packaging away from metal objects. If already unpacked, do not mix or jumble batteries.
19. Remove batteries from the scale if it is not to be used for an extended period of time unless it is for emergency purposes.

SAVE THESE INSTRUCTIONS


WHAT'S IN THE BOX

Elis 1C x1
AAA Battery x3 (preinstalled)
User Manual x1

TECHNICAL SPECIFICATIONS

Product Name: Body Fat Scale
Model: R-A033
Unit: lb / kg / st
Capacity: 0.4-440lb / 0.2-200kg
Division: 0.2lb / 0.05kg
Dimension: 280 x 280 x 25mm / 11 x 11 x 1inch
Battery: 3 x 1.5V AAA (included)

ABOUT YOUR ELIS 1C


Icons

	Calibrating Calibrate the scale before each use.
	Overloaded Maximum capacity: 440lb / 200kg.
	Low Battery Replace the batteries.
	Measuring Body Fat
	Bluetooth Connectivity ✱ Bluetooth is connecting. ✱ Bluetooth is connected.
st lb kg	Units Switched via the app.

* To preserve power, the scale will automatically shut off after 10 seconds of inactivity.

BEFORE USE

Remove the insulation strip from the battery compartment before initial use.

INSTALLING THE APP & PAIRING

1. Download the Renpho Health App

Scan the QR code below or search for “Renpho Health” from the App Store or Google Play to download the Renpho Health App.

Note: Due to continuous updates and improvements, the “Renpho Health” App may appear slightly different.

2. Pair the Scale via Bluetooth

Before pairing, enable Bluetooth on your smartphone.



- iOS devices (version: 12.0 or later) require Bluetooth permission.
- Android devices (version: 7.0 or later) additionally need location permission.

Method 1

- Open the Renpho Health App and step onto the scale's electrodes with bare, dry feet.




- Tap "Connect" to pair when a window pops up on the app.

Method 2


- Open the Renpho Health App and step onto the scale's electrodes with bare, dry feet.
- Tap "+" in the upper-right corner of the "Device" page to connect manually.
- Tap the icon of the Body Fat Scale R-A033 once it pops up.

CALIBRATING YOUR ELIS 1C

Calibrate the scale before each use to ensure accuracy.

1. Place the scale on a flat, hard surface. Avoid any carpeting.
2. Step on the scale with a single foot to activate it and then step off.
3. The display will show "CAL" and then "0.0" or "0.00" to indicate calibrating is complete.


CAL → 00

or

CAL → 0.00


USING YOUR ELIS 1C

1. Make sure the scale is placed on a flat, hard surface.
2. After calibrating the scale, step on the scale's electrodes with bare, dry feet, keep still and wait for the final measurement to be shown on the scale's display.
3. After measuring, the scale will transmit the data to the app via Bluetooth. Open the app to sync and view your detailed measurements, track your progress, and access additional features and insights.

REPLACING THE BATTERIES

Always insert batteries correctly according to the polarity (+ and -) marked on the battery and scale.

Note:

- Clean the battery contacts and also those of the scale prior to battery installation.
- Always select the correct size and type of battery most suitable for the intended use. Do not mix new and used batteries or batteries of different types or brands.
- Replace all batteries of a set at the same time.

TROUBLESHOOTING GUIDE

1. The scale cannot detect body fat and/or other data.

- a. Ensure your feet are dry and bare.
- b. Stand and keep still until the final measurement appears on the scale's display.

If still not available, try applying moisturizer to your feet and measure again or ask another person to try the 2 steps mentioned above.

2. I registered an account on the RENPHO website but failed to use it to log in to the Renpho Health App.

The RENPHO website and Renpho Health App are two separate systems. Please create an account for the Renpho Health App.

WARRANTY POLICY

Your RENPHO product purchase is covered by a one-year limited manufacturer warranty from the date of delivery.

For warranty terms and conditions, please visit:

<https://renpho.com/pages/warranty-terms-and-conditions>

Note: Product registration is not required for the warranty.

If you choose not to register your product, it will not void the product warranty.

CUSTOMER SERVICE

Please feel free to contact us if you have any questions or concerns. RENPHO Customer Service Team guarantees a quick response and hassle-free solutions to any issue you may have within business hours.

✉ **Email: ventas@360mobile.es**

☎ **TELF: +34 910 254 378**

* For defective products or the return of items, please contact us with your order number within the specified warranty period. DO NOT dispose of any product parts as they may be required for inspection/repair.

Notice:

This marking indicates that this product should not be disposed with other household wastes throughout the EU. To prevent possible harm to the environment or human health from uncontrolled waste disposal, recycle it responsibly to promote the sustainable reuse of material resources. To return your used device, please use the return and collection systems or contact the retailer where the product was purchased. They can take this product for environmental safe recycling.

FCC REGULATORY COMPLIANCE POLICY

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

(1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Warning: Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications.

However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

RF Exposure Compliance

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

ISED Regulatory compliance

This device contains licence-exempt transmitters that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s).

Operation is subject to the following two conditions:

(1) This device may not cause interference.

(2) This device must accept any interference, including interference that may cause undesired operation of the device.

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This equipment complies with IC RSS-102 radiation exposure limits set forth for an uncontrolled environment.

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Smart Health Simplified

UK Importer: EEON(UK)LIMITED

Kemp House 160 City Road, London, EC1V 2NX, United Kingdom

EU Importer: Renpho GmbH

Lehmweg 17, 20251 Hamburg Hoheluft-Ost, Deutschland, Germany
Tel: +49 40 23969285 (DE ONLY)

Manufacturer: Shenzhen Ruiyi Business Technology Co., Ltd.

No. 810-C063, 8th Floor, Xiangbin International Financial Centre,
No.18, West Free Trade Street, China Special Economic Zone, Qianhai
Bay, Shenzhen, Guangdong Province, 518000 China

UK REP Richers & Vision Ltd.

Lion Court, Farnham Road, Bordon, England, GU35 0NF

EU REP C&E Connection E-Commerce (DE) GmbH

Zum Linnegraben 20, 65933, Frankfurt am Main, Germany
Email: info@ce-connection.de

*Please have your order number ready before contacting customer support.

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