RENPHO

USER MANUAL

Terea Lite X

Model: R-WMF07

CONTENTS

• ENGLISH	01~10
-----------	-------

• DEUTSCH	11~20
-----------	-------

• FRANÇAIS	21~30
------------	-------

• ITALIANO	31~40
------------	-------

• ESPAÑOL	41~50
-----------	-------

IMPORTANT SAFETY INSTRUCTIONS

READ ALL INSTRUCTIONS BEFORE USING.

1. This device is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the device by a person responsible for their safety.
2. Children should be supervised to ensure that they do not play with the device.
3. **WARNING** – To reduce the risk of burns, fire, electric shock, or injury to persons.
 - a) This device should never be left unattended when plugged in. Unplug from outlet when not in use, and before putting on or taking off parts.
 - b) Do not operate under blanket or pillow. Excessive heating can occur and cause fire, electric shock, or injury to persons.
 - c) Close supervision is necessary when this device is used by, on, or near children, invalids, or disabled persons.
 - d) Use this device only for its intended use as described in this manual. Do not use attachments not recommended by the manufacturer.
 - e) Never operate this device if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water.
 - f) Do not carry this device by supply cord or use cord as a handle.
 - g) Keep the cord away from heated surfaces.
 - h) Never operate the device with the air openings blocked. Keep the air openings free of lint, hair, etc..

- i) Never drop or insert any object into any opening.
- j) Do not use the device outdoors.
- K) Do not operate where aerosol (spray) products are being used or where oxygen is being administered.
- l) To disconnect, turn all controls to the off position, then remove plug from outlet.
- 4. The device is only to be used with the power supply unit provided with the device.
- 5. The device has a heated surface. Persons insensitive to heat must be careful when using the device.
- 6. **WARNING:** Use heated surfaces carefully. May cause serious burns. Do not use over insensitive skin areas or in the presence of poor circulation. The unattended use of heat by children or incapacitated persons may be dangerous.
- 7. **WARNING:** Do not stand on the device. Use only while seated.
- 8. The device is intended for household use only and not intended for medical use. It should not be used in treatment or alleviation of disease or abnormal physical state.

SAVE THESE INSTRUCTIONS

ABOUT THE PRODUCT

1. What's in the Box

Terea Lite X × 1

Power Adapter × 1

2. Technical Specifications

Product Name: Foot Massager

Model: R-WMF07

Input: 12V $\equiv$ 2.5A

Rated Power: 30W

Power Cord Length: 5.9 ft / 1.8 m

Product Weight: 5.4 lb / 2.45 kg


Product Size: 14 × 12.2 × 7 in /
355 × 310 × 177 mm

3. Product Overview

Control Panel


Detachable Foot
Cover with Zipper

Power Port

Note: Do not unzip the foot massager.

4. Control Panel Instruction

USING THE PRODUCT

*Do not use the device for more than 30mins each time.

*Do not stand up while massaging.

1. Connect the device to the power adapter and plug it into a power outlet. Once the device is plugged in, the power indicator will flash white, indicating that it is in standby mode.
2. The device must be used with the power adapter provided.

Buttons Instruction



Button	Operation
	Press to turn ON / OFF the device. Press and hold to enable / disable the heat. [Enable: Solid red (default)]
	Press to cycle through 15 / 30 min. - 15 min: Solid white (Default) - 30 min: Solid red
	Press to cycle through kneading speeds: - Low: Solid white (Default) - High: Solid red - Off: Off Press and hold to cycle through kneading directions: - Outward (Default) - Inward
	Press to cycle through pressure intensities: - Low: Solid white (Default) - Medium: Solid green - High: Solid red - OFF: Off Press and hold to cycle through pressure modes: - Insteps & heels (Default) - Insteps only - Heels only

Note:

1. Unplug the device after each use.
2. This device has three functions: heat, kneading, and compression. Each function can be used separately or in combination.
3. Memory Function: The massager automatically recalls the last selected settings when you press the Power Button again. However, the memory function will not work if the power is cut off.

MAINTENANCE AND CLEANING

1. It is prohibited for non-professional personnel to repair the device or replace its parts.
2. Place and store the device in a cool and dry environment. Do not expose the device to direct sunlight. Keep away from humidity, heat, and/or corrosive gas.
3. Do not place any heavy objects on the device.
4. Keep the device clean. Unplug the device before cleaning.
5. Use a soft, damp cloth to wipe the main body for daily cleaning. To wipe out the stains caused by oil or cosmetics, use a soft, damp cloth with neutral detergent or special detergent for leather. Never immerse the device in water or other liquids.
6. The foot covers are detachable and washable. Unzip and wash them. Do not use the device if the covers are detached.

TROUBLESHOOTING GUIDE

Problems	Causes	Solutions
The device does not work.	The power plug is not connected properly.	Reconnect the plug to the outlet and check the power indicator.
	Not using the supplied adapter, or the interface is wrong, or the power supply is insufficient.	Use the adapter provided with the device.
	Applying excessive pressure could cause nodes to stall.	Do not tread on nodes.

Problems	Causes	Solutions
Insufficient heat or no heat.	The heat function is disabled during massage.	Press and hold the Power Button to enable the heat function.
	The device takes time to heat up gradually to avoid discomfort caused by rapid heat.	Wait several minutes after turning on the device.
	Thick socks are worn during massage, reducing heat sensation.	Wear thinner socks.
The device shuts off suddenly.	The plug is loose, causing poor contact.	Re-plug the adapter and restart the device.
	The timer runs out.	Restart the device for another session.

WARRANTY POLICY

Your RENPHO product purchase is covered by a one-year limited manufacturer warranty from the date of delivery.

For warranty terms and conditions, please visit:

<https://renpho.com/pages/warranty-terms-and-conditions>

Note: Product registration is not required for the warranty. If you choose not to register your product, it will not diminish the product warranty.

CUSTOMER SERVICE

Please feel free to contact us if you have any questions or concerns. RENPHO Customer Service Team guarantees a quick response and hassle-free solutions to any issue you may have within business hours.

TEL: +1(844) 417 0149 (US&CA)

Monday-Friday 9:00AM-4:30PM

1800 MY RENPHO (1800 69 7367) (AU ONLY)

Monday-Friday 9:00AM-6:00PM (AEST)

+44(0)7434 666088 (UK ONLY)

Monday-Friday 9:00AM-5:00PM (GMT)

Email: support@renpho.com (US&CA)

support-au@renpho.com (AU)

support-eu@renpho.com (UK&EU)

*For defective products or the return of items, please contact us with your order number within the specified warranty period. DO NOT dispose of any product parts as they may be required for inspection/repair.

1. Meaning of the crossed-out wheeled bin symbol

The crossed-out wheeled bin symbol means that you are legally obliged to dispose of these appliances separately from unsorted municipal waste, i.e. old appliances must be disposed of separately from household waste. Batteries or rechargeable batteries that are not permanently installed must be removed beforehand and disposed of separately.

2. Options for returning old appliances

Owners of WEEE can return or collect it within the framework of the possibilities set up and available by public waste management authorities to ensure proper disposal.

Hereby, Shenzhen Ruiyi Business Technology Co., Ltd. declares that this equipment complies with Directive 2014/30/EU, Directive 2014/35/EU, Directive 2011/65/EU, The Electrical Equipment (Safety) Regulations 2016, The Electromagnetic Compatibility Regulations 2016, Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012, and other relevant provisions. The full text of the Declaration of Conformity can be obtained at the following internet address: <https://renpho.eu/pages/compliance-data>

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<https://renpho.eu/pages/compliance-data>

FCC Regulatory Compliance

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Warning: changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

ISED Regulatory Compliance

This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s). Operation is subject to the following two conditions:

- (1) This device may not cause interference.
- (2) This device must accept any interference, including interference that may cause undesired operation of the device.

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